

## BME2 – Biomedical Ultrasonics

### *Lecture 8 Bioeffects II- Mechanical Effects*

*Dr. Constantin-C. Coussios*



Office: 43 Banbury Road

Tel: 01865-(2)74750

Email: [constantin.coussios@eng.ox.ac.uk](mailto:constantin.coussios@eng.ox.ac.uk)

*Acknowledgments: Prof. Ronald A. Roy, George Eastman Visiting Professor 2006-07*



Department of Engineering Science

## **Contents**

- 8.1. Acoustic streaming
- 8.2. Radiation force and standing waves
- 8.3. Mechanical effects of inertial and stable cavitation
- 8.6. Case study: lithotripsy
- 8.4. Case study: ultrasound-enhanced thrombolysis

## 8.1 Acoustic streaming [see Hamilton & Blackstock, Chapter 7]

- In a region of fluid where a sound field exists, the pressure and velocity vary with time.
- In general, the temporal averages of these quantities are not zero.
- The average over time of the velocity is called *acoustic streaming*.



- Streaming is strongest
  - At high frequencies
  - In highly attenuative media
  - In highly non-linear media
  - At high pressure amplitudes
  - Under pulsed excitation (for a given acoustic power)
- In water at room temperature, for a diagnostic ultrasound pulse of centre frequency 3.5 MHz, streaming velocities are about **0.33 mm/s** for a 0.1 W acoustic input.
- By maintaining the same acoustic power input (0.1 W) and by pulsing the ultrasound, the streaming velocity can be increased to 10 cm/s at the focal distance!

## 8.2 Acoustic Radiation Force


- A particle placed in an ultrasound field will scatter and absorb part of the incident field, resulting in a net momentum loss over the dimension of the particle along the acoustic axis.

- This results in the particle experiencing a net acoustic radiation force, given by

$$F_A = - \left( \frac{\pi p_0^2 V_c \kappa_1}{2\lambda} \right) \cdot \left[ \frac{5\rho_2 - 2\rho_1}{2\rho_2 + \rho_1} - \frac{\kappa_2}{\kappa_1} \right] \cdot \sin \left( \frac{2\pi}{\lambda} x \right)$$



where  $V_c$  is the particle volume,  $\lambda$  the wavelength,  $p_0$  the pressure amplitude,  $x$  the distance along the acoustic axis, and  $(\rho_2, \kappa_2)$  and  $(\rho_1, \kappa_1)$  are the density and compressibility of the particle and surrounding medium respectively.

- When harnessed within a standing wave field, acoustic forces can be used to separate biological particles based on their density and compressibility contrast!



### 8.3 Mechanical effects of inertial and stable cavitation

- A bubble collapsing inertially near a boundary will cause jetting:



- Such bubble collapses near a cell wall are one of the suspected mechanisms for *sonoporation* (increased permeability of cellular membranes to large molecules).

*(Photos and video courtesy of Lawrence A Crum, Univ. of Washington)*

- A bubble undergoing stable oscillations will generate significant amounts of *microstreaming*





*(Video courtesy of Eleanor Stride, Univ. College London)*

➤ CASE STUDY I:

## EXTRACORPOREAL SHOCK-WAVE LITHOTRIPSY



- Jetting from a cavitation cloud!
- One of several physical mechanisms for kidney stone destruction




*(Photos and video courtesy of Michael Bailey, Univ. of Washington)*

➤ CASE STUDY I: **EXTRACORPOREAL SHOCK-WAVE LITHOTRIPSY**


- Jetting from a cavitation cloud!
- One of several physical mechanisms for kidney stone destruction



*(Photos and video courtesy of Michael Bailey, Univ. of Washington)*



Eroding an artificial kidney stone...



*(Photos and video courtesy of Michael Bailey, Univ. of Washington)*




## ➤ CASE STUDY II:

### ULTRASOUND-ENHANCED THROMBOLYSIS

#### Ultrasound-Enhanced Thrombolysis: Clinical Motivation

- Blood clots and thrombi are largely responsible for ischaemic stroke thrombotic events.
- Thrombolytic action of **rt-PA**




- Systemic delivery of rt-PA can cause internal bleeding.
- **Localized activation key to successful & safe clinical use.**

# Thrombolysis Study: Clot Formation & Exposure

- Whole Blood Clots:
  - 2-3 ml aliquots of whole pig blood clot in test tube
  - Incubated at 37 C for 3 hours
  - Refrigerated 3 days prior to use
- Pulsed or CW Ultrasound Exposures
  - 120 KHz center frequency
  - PRP = 600  $\mu$ s
  - Exposure Duration = 30 min.
  - Duty Cycle from 0% to CW
- Clot added to rt-PA + Plasma
- Weighed before & after exposure




## Apparatus for Cavitation Detection







# Excitation regimes – cavitation study



## Correlation of cavitation with thrombolysis



## Observations of clot surface



Labeled with mouse anti-fibrinogen antibody and stained with horseradish peroxidase-linked goat anti-mouse IgG

## Conclusions: the future of biomedical ultrasound

- Ultrasound is cheap, portable, presents no risk to the user and little or no risk to the patient.
- Ultrasound imaging has become widely accepted in the clinical arena and is especially well-suited to imaging of high acoustic contrasts (suffers by comparison to MRI for image resolution).
- The future of ultrasound applied to therapy is bright: there are numerous emerging and under-explored applications, requiring deeper understanding of the interactions of sound (and bubbles) with tissue and cells